

## Product Datasheet

### Human GRIA2 ELISA Kit (orb1199002)

<b>Description</b>	Human Glutamate Receptor, Ionotropic, AMPA 2 (GRIA2) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.312-20ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">P42262</a>
<b>Sensitivity</b>	0.115ng/mL
<b>Entrez</b>	<a href="#">2891</a>
<b>Expiration Date</b>	Please enquire.