

Product Datasheet

Human ELOVL1 ELISA Kit (orb1199138)

Description	Human Elongation Of Very Long Chain Fatty Acids Like Protein 1 (ELOVL1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9BW60
Sensitivity	0.056ng/mL
Entrez	64834
Expiration Date	Please enquire.