

Product Datasheet

Human CRHR1 ELISA Kit (orb1199288)

Description	Human Corticotropin Releasing Hormone Receptor 1 (CRHR1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P34998
Sensitivity	0.060ng/mL
Entrez	1394
Expiration Date	Please enquire.