

Product Datasheet

Human CRELD1 ELISA Kit (orb1199290)

Description	Human Cysteine Rich With EGF Like Domains Protein 1 (CRELD1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q96HD1
Sensitivity	0.113ng/mL
Entrez	78987
Expiration Date	Please enquire.