

## Product Datasheet

### Human CDNF ELISA Kit (orb1199381)

**Description** Human Cerebral Dopamine Neurotrophic Factor (CDNF) ELISA Kit

**Reactivity** Human

**Range** 0.156-10ng/mL

**Tested Applications** ELISA

**Note** For research use only

**Application notes** sandwich

**Sample Types** serum, plasma, tissue homogenates, cell lysates or other biological fluids.

**Uniprot ID** [Q49AH0](#)

**Sensitivity** 0.057ng/mL

**Entrez** [441549](#)

**Expiration Date** Please enquire.