

Product Datasheet

Human CALPRO ELISA Kit (orb1199414)

Description	Human Calprotectin (CALPRO) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P06702, P05109
Sensitivity	13.12pg/mL
Entrez	6280, 6279
Expiration Date	Please enquire.