

Product Datasheet

Human ARNTL ELISA Kit (orb1199491)

Description	Human Aryl Hydrocarbon Receptor Nuclear Translocator Like Protein (ARNTL) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O00327
Sensitivity	0.057ng/mL
Entrez	406
Expiration Date	Please enquire.