

Product Datasheet

Human ADM2 ELISA Kit (orb1199563)

Description	Human Adrenomedullin 2 (ADM2) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q7Z4H4
Sensitivity	26.2pg/mL
Entrez	79924
Expiration Date	Please enquire.