

Product Datasheet

Human ABCC8 ELISA Kit (orb1199584)

Description	Human ATP Binding Cassette Transporter C8 (ABCC8) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q09428
Sensitivity	0.128ng/mL
Entrez	6833
Expiration Date	Please enquire.