

Product Datasheet

Mouse IL1F9 ELISA Kit (orb1199611)

Description	Mouse Interleukin 1 Family, Member 9 (IL1F9) ELISA Kit
Reactivity	Mouse
Range	7.81-500pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q3U0P4
Sensitivity	2.6pg/mL
Entrez	215257
Expiration Date	Please enquire.