

Product Datasheet

Rat TLR4 ELISA Kit (orb1199665)

Description	Rat Toll Like Receptor 4 (TLR4) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9QX05
Sensitivity	0.060ng/mL
Entrez	29260
Expiration Date	Please enquire.