

Product Datasheet

Mouse Pgp ELISA Kit (orb1199787)

Description	Mouse Permeability Glycoprotein (Pgp) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P06795
Sensitivity	0.051ng/mL
Entrez	18669, 5243
Expiration Date	Please enquire.