

Product Datasheet

Rabbit NT-ProBNP ELISA Kit (orb1199818)

Description	Rabbit N-Terminal Pro Brain Natriuretic Peptide (NT-ProBNP) ELISA Kit
Reactivity	Rabbit
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Sensitivity	11.3pg/mL
Expiration Date	Please enquire.