

Product Datasheet

Rat NT3 ELISA Kit (orb1199824)

Description	Rat Neurotrophin 3 (NT3) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P18280
Sensitivity	0.063ng/mL
Entrez	81737
Expiration Date	Please enquire.