

## Product Datasheet

### Rat NPS ELISA Kit (orb1199832)

<b>Description</b>	Rat Neuropeptide S (NPS) ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	15.6-1000pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	competitive inhibition
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<a href="#">POCOP7</a>
<b>Sensitivity</b>	4.75pg/mL
<b>Entrez</b>	<a href="#">100360071</a>
<b>Expiration Date</b>	Please enquire.