

Product Datasheet

Mouse ITGa4 ELISA Kit (orb1199942)

Description	Mouse Integrin Alpha 4 (ITGa4) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q00651
Sensitivity	0.071ng/mL
Entrez	16401
Expiration Date	Please enquire.