

## Product Datasheet

### Mouse GnRHR ELISA Kit (orb1200001)

<b>Description</b>	Mouse Gonadotropin Releasing Hormone Receptor (GnRHR) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.39-25ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q01776</a>
<b>Sensitivity</b>	0.116ng/mL
<b>Entrez</b>	<a href="#">14715</a>
<b>Expiration Date</b>	Please enquire.