

Product Datasheet

Rat GnRH ELISA Kit (orb1200002)

Description	Rat Gonadotropin Releasing Hormone (GnRH) ELISA Kit
Reactivity	Rat
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P07490
Sensitivity	29pg/mL
Entrez	25194
Expiration Date	Please enquire.