

## Product Datasheet

### Rat bEP ELISA Kit (orb1200184)

<b>Description</b>	Rat Beta-Endorphin (bEP) ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	15.6-1000pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	competitive inhibition
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cerebrospinal fluid, cell culture supernates or other biological fluids.
<b>Sensitivity</b>	4.35pg/mL
<b>Expiration Date</b>	Please enquire.