

Product Datasheet

Mouse AGRP ELISA Kit (orb1200236)

Description	Mouse Agouti Related Protein (AGRP) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P56473
Sensitivity	28pg/mL
Entrez	11604
Expiration Date	Please enquire.