

## Product Datasheet

### Rat ADH1 ELISA Kit (orb1200240)

<b>Description</b>	Rat Alcohol Dehydrogenase 1 (ADH1) ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	1.56-100ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<b>P06757</b>
<b>Sensitivity</b>	0.59ng/mL
<b>Entrez</b>	<b>24172</b>
<b>Expiration Date</b>	Please enquire.