

Product Datasheet

Mouse a1AT ELISA Kit (orb1200258)

Description	Mouse Alpha-1-Antitrypsin (a1AT) ELISA Kit
Reactivity	Mouse
Range	3.12-200ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P22599
Sensitivity	1.42ng/mL
Entrez	20701
Expiration Date	Please enquire.