

Product Datasheet

Human IL1F9 ELISA Kit (orb1200294)

Description	Human Interleukin 1 Family, Member 9 (IL1F9) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9NZH8
Sensitivity	5.9pg/mL
Entrez	56300
Expiration Date	Please enquire.