

Product Datasheet

Human IGF2BP3 ELISA Kit (orb1200295)

Description	Human Insulin Like Growth Factor 2 mRNA Binding Protein 3 (IGF2BP3) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O00425
Sensitivity	0.114ng/mL
Entrez	10643
Expiration Date	Please enquire.