

Product Datasheet

Human TLR1 ELISA Kit (orb1200331)

| | |
|----------------------------|---|
| Description | Human Toll Like Receptor 1 (TLR1) ELISA Kit |
| Reactivity | Human |
| Range | 0.156-10ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | tissue homogenates, cell lysates or other biological fluids. |
| Uniprot ID | Q15399 |
| Sensitivity | 0.072ng/mL |
| Entrez | 7096 |
| Expiration Date | Please enquire. |