

Product Datasheet

Human S100A9 ELISA Kit (orb1200377)

Description	Human S100 Calcium Binding Protein A9 (S100A9) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P06702
Sensitivity	0.24ng/mL
Entrez	6280
Expiration Date	Please enquire.