

Product Datasheet

Human IκBa ELISA Kit (orb1200546)

Description	Human Inhibitory Subunit Of NF Kappa B Alpha (IκBa) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P25963
Sensitivity	0.055ng/mL
Entrez	4792
Expiration Date	Please enquire.