

Product Datasheet

Human GHITM ELISA Kit (orb1200588)

Description	Human Growth Hormone Inducible Transmembrane Protein (GHITM) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9H3K2
Sensitivity	0.13ng/mL
Entrez	27069
Expiration Date	Please enquire.