

## Product Datasheet

### Mouse Dyn ELISA Kit (orb1200639)

<b>Description</b>	Mouse Big Dynorphin (Dyn) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	4.68-300pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">O35417</a>
<b>Sensitivity</b>	2.4pg/mL
<b>Entrez</b>	<a href="#">18610</a>
<b>Expiration Date</b>	Please enquire.