

## Product Datasheet

### Human DEFb1 ELISA Kit (orb1200647)

<b>Description</b>	Human Defensin Beta 1 (DEFb1) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	62.5-4000pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<b>P60022</b>
<b>Sensitivity</b>	25.3pg/mL
<b>Entrez</b>	<b>1672</b>
<b>Expiration Date</b>	Please enquire.