

## Product Datasheet

### Human Btk ELISA Kit (orb1200726)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Human Bruton'S Tyrosine Kinase (Btk)<br>ELISA Kit               |
| <b>Reactivity</b>          | Human   |
| <b>Range</b>               | 0.156-10ng/mL   |
| <b>Tested Applications</b> | ELISA   |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | sandwich  |
| <b>Sample Types</b>        | tissue homogenates, cell lysates or<br>other biological fluids. |
| <b>Uniprot ID</b>          | <a href="#">Q06187</a>  |
| <b>Sensitivity</b>         | 0.059ng/mL  |
| <b>Entrez</b>              | <a href="#">695</a>   |
| <b>Expiration Date</b>     | Please enquire.   |