

Product Datasheet

Human a1ACT ELISA Kit (orb1200805)

| | |
|----------------------------|--|
| Description | Human Alpha-1-Antichymotrypsin (a1ACT) ELISA Kit |
| Reactivity | Human |
| Range | 3.12-200ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | serum, plasma or other biological fluids. |
| Uniprot ID | P01011 |
| Sensitivity | 1.42ng/mL |
| Entrez | 12 |
| Expiration Date | Please enquire. |