

## Product Datasheet

### Rat EDN2 ELISA Kit (orb1200827)

<b>Description</b>	Rat Endothelin 2 (EDN2) ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	7.81-500pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">P23943</a>
<b>Sensitivity</b>	2.96pg/mL
<b>Entrez</b>	<a href="#">24324</a>
<b>Expiration Date</b>	Please enquire.