

Product Datasheet

Porcine TGFb1 ELISA Kit (orb1200838)

Description	Porcine Transforming Growth Factor Beta 1 (TGFb1) ELISA Kit
Reactivity	Porcine
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, platelet-poor plasma, tissue homogenates, cell culture supernates or other biological fluids
Uniprot ID	P07200
Sensitivity	10.5pg/mL
Entrez	397078
Expiration Date	Please enquire.