

Product Datasheet

Bovine IP10 ELISA Kit (orb1200843)

| | |
|----------------------------|--|
| Description | Bovine Interferon Gamma Induced Protein 10kDa (IP10) ELISA Kit |
| Reactivity | Bovine |
| Range | 0.156-10ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Uniprot ID | Q2KIQ8 |
| Sensitivity | 0.058ng/mL |
| Entrez | 615107 |
| Expiration Date | Please enquire. |