

## Product Datasheet

### Bovine IL1RA ELISA Kit (orb1200853)

<b>Description</b>	Bovine Interleukin 1 Receptor Antagonist (IL1RA) ELISA Kit
<b>Reactivity</b>	Bovine
<b>Range</b>	31.2-2000pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">O77482</a>
<b>Sensitivity</b>	13.5pg/mL
<b>Entrez</b>	<a href="#">281860</a>
<b>Expiration Date</b>	Please enquire.