

Product Datasheet

Porcine IL1b ELISA Kit (orb1200854)

Description	Porcine Interleukin 1 Beta (IL1b) ELISA Kit
Reactivity	Porcine
Range	62.5-4000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q29082
Sensitivity	26.2pg/mL
Entrez	396565
Expiration Date	Please enquire.