

Product Datasheet

Bovine IL1b ELISA Kit (orb1200855)

Description	Bovine Interleukin 1 Beta (IL1b) ELISA Kit
Reactivity	Bovine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P09428
Sensitivity	6.1pg/mL
Entrez	281251
Expiration Date	Please enquire.