

Product Datasheet

Porcine IL18 ELISA Kit (orb1200858)

Description	Porcine Interleukin 18 (IL18) ELISA Kit
Reactivity	Porcine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O19073
Sensitivity	5.8pg/mL
Entrez	397057
Expiration Date	Please enquire.