

Product Datasheet

Porcine IL17 ELISA Kit (orb1200860)

Description	Porcine Interleukin 17 (IL17) ELISA Kit
Reactivity	Porcine
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q60I29
Sensitivity	12.5pg/mL
Entrez	449530
Expiration Date	Please enquire.