

Product Datasheet

Porcine IL13 ELISA Kit (orb1200866)

Description	Porcine Interleukin 13 (IL13) ELISA Kit
Reactivity	Porcine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q95J68
Sensitivity	5.8pg/mL
Entrez	396721
Expiration Date	Please enquire.