

Product Datasheet

Bovine IFN α ELISA Kit (orb1200889)

Description	Bovine Interferon Alpha (IFN α) ELISA Kit
Reactivity	Bovine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P05007
Sensitivity	5.9pg/mL
Entrez	282257
Expiration Date	Please enquire.