

Product Datasheet

Porcine FGF7 ELISA Kit (orb1200893)

Description	Porcine Fibroblast Growth Factor 7 (FGF7) ELISA Kit
Reactivity	Porcine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9N198
Sensitivity	5.9pg/mL
Entrez	397281
Expiration Date	Please enquire.