

Product Datasheet

Porcine FGF1 ELISA Kit (orb1200897)

Description	Porcine Fibroblast Growth Factor 1, Acidic (FGF1) ELISA Kit
Reactivity	Porcine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	F1RM38
Sensitivity	6.6pg/mL
Entrez	397497
Expiration Date	Please enquire.