

## Product Datasheet

### Porcine BMP2 ELISA Kit (orb1200907)

<b>Description</b>	Porcine Bone Morphogenetic Protein 2 (BMP2) ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	31.2-2000pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">E1AXT2</a>
<b>Sensitivity</b>	11.8pg/mL
<b>Entrez</b>	<a href="#">100157103</a> , <a href="#">100152318</a>
<b>Expiration Date</b>	Please enquire.