

Product Datasheet

Mouse SYP ELISA Kit (orb1200964)

Description	Mouse Synaptophysin (SYP) ELISA Kit
Reactivity	Mouse
Range	0.625-40ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cerebrospinal fluid or other biological fluids.
Uniprot ID	Q62277
Sensitivity	0.278ng/mL
Entrez	20977
Expiration Date	Please enquire.