

Product Datasheet

Mouse NT-ProANP ELISA Kit (orb1201033)

Description	Mouse N-Terminal Pro Atrial Natriuretic Peptide (NT-ProANP) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sensitivity	29pg/mL
Expiration Date	Please enquire.