

## Product Datasheet

### Mouse MIP3b ELISA Kit (orb1201066)

<b>Description</b>	Mouse Macrophage Inflammatory Protein 3 Beta (MIP3b) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	15.6-1000pg/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<a href="#">O70460</a>
<b>Sensitivity</b>	5.9pg/mL
<b>Entrez</b>	<a href="#">24047</a>
<b>Expiration Date</b>	Please enquire.