

Product Datasheet

Mouse MCP2 ELISA Kit (orb1201071)

Description	Mouse Monocyte Chemotactic Protein 2 (MCP2) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9Z121
Sensitivity	6.3pg/mL
Entrez	20307
Expiration Date	Please enquire.