

Product Datasheet

Mouse bTG ELISA Kit (orb1201222)

Description	Mouse Beta-Thromboglobulin (bTG) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P62325
Sensitivity	34pg/mL
Entrez	12226
Expiration Date	Please enquire.